

From Idea to Image™



6.5EC128-EP

2D Endocavity Probe

***Providing high performance
and increased durability...***

Endocavity probes undergo frequent disinfection procedures and are prone to mechanical and chemical damages that can lead to expensive repairs and replacements. To address this issue, we have enhanced the durability of our endocavity probes by using our rigid Encap™ technology lens.

Features

- 6.5MHz center frequency
- High performance Encap™ acoustic design
- 128 element micro-convex array
- 10mm radius of curvature
- Ergonomic handle design
- Disposable needle guide

Sound Technology, Inc.

1363 South Atherton Street
State College, PA 16803
USA

T: 814.234.4377
F: 814.234.5033
info@sti-ultrasound.com

a subsidiary of

ANALOGIC

For information on how **STI's** products and solutions can improve your image in the ultrasound market, please contact us at info@sti-ultrasound.com or visit our website at sti-ultrasound.com.